


SlideShaft®

Adjustable Shaft Length Neuro-stimulation Probe

- All Neuro-stimulating Probe Shaft Sizes in One Probe (130mm to 270mm)
- Adjust Probe Shaft Length to Surgical Procedure
- Realtime Intraoperative Customization of Probe Shaft Length According to Patient Anatomy



Ordering Information: SlideShaft® Stimulator Probes

Item Code	Description
NP001	SlideShaft®, Single Use Disposable Flush Tip Probe; 0.2mm tip; 2.2mm shaft with adjustable shaft length 130mm to 270mm electrode; 1.9m wire with DIN 42 802 connector
NP002	SlideShaft®, Single Use Disposable Ball Tip Probe; 2.3mm ball with adjustable shaft length 130mm to 270mm electrode; 1.9m wire with DIN 42 802 connector



NeuralPathways, LLC

NeuralPathways, LLC
2674 East Main Street, #432, Ventura, CA, 93003
Phone: 805.290.2710
www.neuralpathways.com
e-mail: info@neuralpathways.com

Manufactured by:
Technomed Europe
Amerikalaan 71
6199AE Maastricht Airport
The Netherlands

Manufactured Exclusively for:
NeuralPathways, LLC
2674 East Main Street
Ventura, CA 93003 USA