


Your Australian Supplier for Concentric EMG Needle Electrodes and Cables



Meditrend
INTERNATIONAL

- High quality concentric EMG needle electrode made for the European Market now available in Australia
- Platinum core recording electrode with ultra-sharp Trocar point ground at a unique 12° Angle.
- Requires reusable cable (F Series or M Series @ \$50).
- Sold in packages of twenty five (25) needles.
- Pre-sterilized. Single use only.
- Enquire about any other Needle gauge sizes

+61.411 519 104 (Trevor)
+61.412 298 065 (Melville)

sales@meditrend.com.au

Visit us at
www.meditrend.com.au

From
\$ 159
+ GST
For 25 Needles

FREE Postage in
Australia

Silver Series (Spec Nr. 15.1)

Stainless Wire	\$ 169.00 + GST
Silver Wire	\$ 179.00 + GST
Platinum Wire	\$ 199.00 + GST

Plastic Series (Spec Nr. 10.1)

Stainless Wire	\$ 159.00 + GST
Silver Wire	\$ 179.00 + GST
Platinum Wire	\$ 189.00 + GST